

Epithelium 26[®] LINGERIE "CUSHIONS"



BACKGROUND:

What is the ball of the foot?

The foot has a natural weight distributor called the "plantar ball". This sub-cutaneous tissue is 1 cm thick and made up of many fatty vesicles that deform under pressure. Unfortunately, with age there is often an irreversible change in the plantar ball which causes intense pain at the front of the foot, as well as the formation of calluses.

INDICATIONS :

- Pressure on the ball of the foot
- Pain on the ball of the foot, calluses, etc.
- Rheumatism, algodystrophy, Morton disease
- Prolonged standing, wearing high heels
- Intense activity or sport involving the forefoot

PROPERTIES :

- The Epithelium 26[®] "lingerie" cushion is 2.2mm thick and acts as a substitute plantar sole. The Epithelium 26[®] module located under the forefoot distributes pressure in painful areas.
- It combines the effectiveness of traditional plantar cushions with the aesthetic appeal of lingerie due to its white fabric/material with embroidery effect.
- The fabric was selected for its ability to evacuate perspiration in addition to its softness, making this cushion very pleasant to wear.

RECOMMENDATIONS FOR USE:

- Wear directly on the skin. Wear socks, tights or stockings over the device for optimal comfort.
- Machine-washable at 40°C or by hand without rubbing.
- Average duration of use: several months



References	Size	Shoe size			Unit
		EU	UK	US	
CC2621	Small	36-38	3.5-5	W5½-7	One pair
CC2622	Medium	39-41	5.5-7.5	W7½-9½ / M5½ -9½	One pair
CC2623	Large	42-45	8-10.5	W10-12½ / M8-10½	One pair