

PLANTAR CUSHIONS with Epithelium 26®



BACKGROUND:

What is the ball of the foot?

The foot has a natural weight distributor called the “plantar ball”. This sub-cutaneous tissue is 1 cm thick and made up of many fatty vesicles that deform under pressure. Unfortunately, with age, there is often an irreversible change in the plantar ball which causes intense pain at the front of the foot, as well as the formation of calluses.

INDICATIONS :

- Pressure on the ball of the foot
- Pain on the ball of the foot, calluses, etc.
- Rheumatism, algodystrophy, Morton disease
- Prolonged standing, wearing high heels
- Intense activity or sport involving the forefoot

PROPERTIES :

- The Epithelium 26® cushion is 2.2mm thick and acts as a substitute plantar sole. The Epithelium 26® module located under the forefoot distributes pressure in painful areas.
- These cushions can be discreetly worn in your regular shoes.
- The fabric was selected for its ability to evacuate perspiration in addition to its softness, making this cushion very pleasant to wear.



RECOMMENDATIONS FOR USE:

- Wear directly against the skin. Wear socks, tights or stockings over the device for optimal comfort.
- Can be machine-washed at 40°C in the washing bag.
- Average duration of use: several months.



Supplied with washing bag



References	Size	Shoe size			Unit
		EU	UK	US	
CC 2612	Small	36-38	3.5-5	W5½-7	One pair
CC 2613	Medium	39-41	5.5-7.5	W7½-9½ / M5½ -7½	One pair
CC 2614	Large	42-45	8-10.5	W10-12½ / M8-10½	One pair