



SUMMER PLANTAR CUSHIONS With Epithelium 26®



BACKGROUND:

What is the ball of the foot?

The foot has a natural weight distributor called the “plantar ball”. This sub-cutaneous tissue is 1 cm thick and made up of many fatty vesicles that deform under pressure. Unfortunately, with age there is often an irreversible change in the plantar ball which causes intense pain at the front of the foot, as well as the formation of calluses.

INDICATIONS:

- Pressure on the ball of the foot
- Pain on the ball of the foot, calluses, etc.
- Rheumatism, algodystrophy, Morton disease
- Prolonged standing, wearing high heels
- Intense activity or sport involving the forefoot

PROPERTIES:

- The Epithelium 26® module located under the forefoot distributes pressure in painful areas.
- Discreet and thin design. Specially designed to be worn with open-toed shoes or heels.
- Summer shoe cushions are maintained in place in the shoe with an adhesive.

RECOMMENDATIONS FOR USE:

- Wear directly against the skin.
- Slide the cushions to the front of the shoe, remove the adhesive seal and attach them directly to the shoe.
- Machine-washable at 40°C or by hand without rubbing.
- Average duration of use: several months.



References	Shoe size			Unit
	EU	UK	US	
0671	34-35	2-2.5	W4-4½	One pair
0672	36-37	3.5-4.5	W5½-6½	One pair
0673	38-39	5-5.5	W7-8 / M5½-6	One pair
0674	40-41	6.5-7.5	W8½-9½ / M6½-7½	One pair
0675	42-43	8-9	W10-11 / M8-9	One pair
0676	44-45	9.5-10.5	W11½-12½ / M9½-10½	One pair