



DIGITUBES®
with Epithelium 26®



Epithelium 26® on half of a Digitube relieves pressure and rubbing.

INDICATIONS:

- Protection of corns and interdigital corns after treatment.

PROPERTIES:

- 1 mm of Epithelium 26® on half of the Digitube® is enough to relieve pressure and prevent the formation of **dorsal and interdigital/soft corns**.
- The chosen material's softness and elasticity makes the device comfortable to wear.
- Their small size makes them easy to wear on neighboring toes.
- Easy to put on and maintain in place.



RECOMMENDATIONS FOR USE:

- Hand-washable
- With only 1.5 cm of Digitube®, a corn is protected for 10 days on average.



References	Size	Unit	Unit
0221	Small	22 mm / 0.88"	Pack of 2 X 20 cm / 8"
0222	Medium	25 mm / 1"	Pack of 2 X 20 cm / 8"
0223	Large	33 mm / 1.32"	Pack of 2 X 20 cm / 8"
0224	S+M+L	22+25+33 mm / 0.88"+1"+1.32"	Pack of 3 X 20 cm / 8"

Consumer Products (with Sales Packaging)

DI261	Small	22 mm / 0.88"	Pack of 3 X10cm/4"
DI262	Medium	25 mm / 1"	Pack of 3x10cl/4"
DI263	Large	33 mm / 1.32"	Pack of 3x10cm/4"
DI264	S+M+L	22+25+33 mm / 0.88"+1"+1.32"	Pack of 3x10cm/4"