



## TOE TIP PROTECTIONS FOR HAMMER TOES with Epithelium 26®

### BACKGROUND :

What is a “hammer toe”?

**A hammer toe is a lifted toe. This problem mainly affects the 2nd, 3rd and 4th toes. This condition is often associated with a “pes cavus”: a foot with an excessively high arch. Folded toes are then crowded between the tops of shoes and the inner soles. Pressure from shoes can create corns on the tops of toe joints (dorsal corn) and on the tips of toes (soft corn).**

### INDICATIONS :

- Hammer toes
- Claw toes
- Dorsal corns

### PROPERTIES :

- When the Epithelium 26® toe protector is placed on top of toes, it distributes pressure and limits the formation of corns on the tops of toes.
- These protectors protect all of the toes, yet are small and easy to wear: the integrated protective gel is only 1 mm thick.
- The fabric is soft and absorbs perspiration, making this device very comfortable to wear. It also provides warmth to the tips of toes, even with poor circulation.



### RECOMMENDATIONS FOR USE:

- Can be machine-washed at 40°C.
- Average duration of use: several months.



washing bag supplied



References	Size	Shoe size			Unit
		EU	UK	US	
0451	Small	36-38	3.5-5	W5½-7	One pair
0452	Medium	39-41	5.5-7.5	W7½-9½ / M5½ -7½	One pair
0453	Large	42-45	8-10.5	W10-12½ / M8-10½	One pair